BRUGSINFORMATIONER DK

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## CARE FLOW

### Urin-Drainage-Systeme - unsteril -

### Urinary-Drainage-Systems - non-sterile -

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INSTRUCTIONS FOR USE

**General Information:**  
Dear user, please read all of the following information carefully before use. Please refer to the table for the specific product features of the individual systems.

**Intended Purpose/Description:**  
Urine drainage systems and bags for connection to catheters and urinary sheaths. Non-sterile urine drainage systems and bags labelled US/UF are only suitable for use with urinary sheaths and for connection to leg bags.

| Product                     | filling volume in ml | Tube length in cm | Leg bag            | Drainage tap/valve | Non-woven coating | Fastening and combination options |                    | <div>CE</div> <div>1<br/>Items suitable for combination are not included in the scope of delivery. These are available to order separately from us.</div> |
|-----------------------------|----------------------|-------------------|--------------------|--------------------|-------------------|-----------------------------------|--------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>CareFlow</b>             |                      |                   | Number of chambers |                    | Reverse side      | Twin holder                       | Elastic band/ Cuff |                                                                                                                                                           |
| UB 2000 US                  | 2000                 | 90/130            |                    | Tap                |                   | x                                 |                    |                                                                                                                                                           |
| UB 2000 USo/<br>SB 2000 USo | 2000                 | 90                |                    |                    |                   | x                                 |                    |                                                                                                                                                           |
| SB 1500 USo/<br>UB 1500 USo | 1500                 | 10/90             |                    |                    |                   | x                                 |                    |                                                                                                                                                           |
| UB 750 US                   | 750                  | 10/30/50          | 1                  | Tap                | x                 | x                                 | x                  |                                                                                                                                                           |
| UB 750 UF                   | 750                  | 50                | 1                  | Valve              | x                 | x                                 | x                  |                                                                                                                                                           |
| UB 500 US                   | 500                  | 10/30/50          | 3                  | Tap                | x                 | x                                 | x                  |                                                                                                                                                           |

- Indication/Application Area:**  
Intended for use on patients with a recumbent indwelling bladder catheter, percutaneous renal fistula or urinary sheath for collecting the drained urine by means of a urine drainage system.
- The patient must be cared for in a way that is tailored to their individual needs. If the attending physician has given any previous instructions, then these must be observed.
- Contraindication and Restriction of Use:**  
We are not aware of any absolute contraindications. Restrictions of use may result from the clinical picture of the patient. Leg bag systems should not be used on patients with macrohematuria. Our urine drainage systems and their tubes must not come into contact with areas of the skin presenting pathological changes.
- Application/Execution/Hygiene:**  
Before manipulating the urine drainage systems in any way (connection or change of the system/urine drainage), please carry out hygienic hand disinfection (wash and disinfect hands). Please wear disposable gloves.
- Take the bag out of the packaging
  - Remove protective cap at cone. If necessary, close non-return valve at drainage tap (if open).
  - When using urinary sheaths, the urine bag must also be discarded during the daily change.
  - Connect the bag system with the urinary sheath by inserting the cone into the funnel or into the silicone connector of the leg bag.
  - Achten Sie darauf, dass die Verpackungseinheit ungeöffnet und unbeschädigt ist und nicht länger als auf dem Etikett angegeben eingesetzt wird.
  - Leg bag: For optimal drainage, ensure the bag is placed below the bladder.
- Warnings/Safety Precautions:**  
Please note that urine drainage systems connected to an indwelling bladder catheter must not be pulled in any way in order to avoid disturbing the position of the catheter and impairing its drainage performance (RISK OF INJURY!) If the bag or medical device is defective and/or no urine enters the bag, discard it immediately and use a new one. If patients affected by conditions that require urinary diversion are fitted with our urine drainage systems, appropriate care/supervision must be ensured. The same applies to small children. Do not reuse, reprocess. Reuse, reprocessing may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in injury, illness or death of the patient. Reuse, reprocessing may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. For device complaints, please contact your supplier. Please also report any serious incident that has occurred in relation to the device to the manufacturer and the competent authority at the user's place of residence.
- Exclusion of Liability:**  
In general, warranty only applies if the connection of approved medical products from GHC to the product you have in hand is carried out correctly and in accordance with the instructions for use. If you prefer to connect these to products from other manufacturers, this is possible, but GHC assumes no liability and any warranty claims become void.

GHC assumes no liability whatsoever for injury or other damage resulting from use of this product other than in accordance with current GHC recommendations. PLEASE NOTE; All information has been issued to the best of our knowledge at the time of printing. Subject to technical changes.

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