BRUGSINFORMATIONER DK

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## CARE FLOW

### Urin-Drainage-Systeme - unsteril -

#### Urinary-Drainage-Systems - non-sterile -

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INSTRUCTIONS FOR USE

**General Information:**  
Dear user, please read all of the following information carefully before use. Please refer to the table for the specific product features of the individual systems.

**Intended Purpose/Description:**  
Urine drainage systems and bags for connection to catheters and urinary sheaths. Non-sterile urine drainage systems and bags labelled US/UF are only suitable for use with urinary sheaths and for connection to leg bags.

| Product                     | filling volume<br>in ml | Tube length<br>in cm | Leg bag               | Drainage<br>tap/valve | Non-woven<br>coating | Fastening and<br>combination options |                       | <br><br>1<br>Items<br>suitable for<br>combination<br>are not<br>included in<br>the scope of<br>delivery. These<br>are available<br>to order<br>separately<br>from us. |
|-----------------------------|-------------------------|----------------------|-----------------------|-----------------------|----------------------|--------------------------------------|-----------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| CareFlow                    |                         |                      | Number of<br>chambers |                       | Reverse side         | Twin holder                          | Elastic band/<br>Cuff |                                                                                                                                                                                                                                                          |
| UB 2000 US                  | 2000                    | 90/130               |                       | Tap                   |                      | x                                    |                       |                                                                                                                                                                                                                                                          |
| UB 2000 USo/<br>SB 2000 USo | 2000                    | 90                   |                       |                       |                      | x                                    |                       |                                                                                                                                                                                                                                                          |
| SB 1500 USo/<br>UB 1500 USo | 1500                    | 10/90                |                       |                       |                      | x                                    |                       |                                                                                                                                                                                                                                                          |
| UB 750 US                   | 750                     | 10/30/50             | 1                     | Tap                   | x                    | x                                    | x                     |                                                                                                                                                                                                                                                          |
| UB 750 UF                   | 750                     | 50                   | 1                     | Valve                 | x                    | x                                    | x                     |                                                                                                                                                                                                                                                          |
| UB 500 US                   | 500                     | 10/30/50             | 3                     | Tap                   | x                    | x                                    | x                     |                                                                                                                                                                                                                                                          |

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- Warnings/Safety Precautions:**  
Please note that urine drainage systems connected to an indwelling bladder catheter must not be pulled in any way in order to avoid disturbing the position of the catheter and impairing its drainage performance (RISK OF INJURY!) If the bag or medical device is defective and/or no urine enters the bag, discard it immediately and use a new one. If patients affected by conditions that require urinary diversion are fitted with our urine drainage systems, appropriate care/supervision must be ensured. The same applies to small children. Do not reuse, reprocess. Reuse, reprocessing may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in injury, illness or death of the patient. Reuse, reprocessing may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. For device complaints, please contact your supplier. Please also report any serious incident that has occurred in relation to the device to the manufacturer and the competent authority at the user's place of residence.
- Exclusion of Liability:**  
In general, warranty only applies if the connection of approved medical products from GHC to the product you have in hand is carried out correctly and in accordance with the instructions for use. If you prefer to connect these to products from other manufacturers, this is possible, but GHC assumes no liability and any warranty claims become void.

GHC assumes no liability whatsoever for injury or other damage resulting from use of this product other than in accordance with current GHC recommendations. PLEASE NOTE; All information has been issued to the best of our knowledge at the time of printing. Subject to technical changes.

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