BRUGSINFORMATIONER DK

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CARE FLOW

Urin-Drainage-Systeme - unsteril -

Urinary-Drainage-Systems - non-sterile -

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INSTRUCTIONS FOR USE

General Information:
Dear user, please read all of the following information carefully before use. Please refer to the table for the specific product features of the individual systems.

Intended Purpose/Description:
Urine drainage systems and bags for connection to catheters and urinary sheaths. Non-sterile urine drainage systems and bags labelled US/UF are only suitable for use with urinary sheaths and for connection to leg bags.

Product	filling volume in ml	Tube length in cm	Leg bag	Drainage tap/valve	Non-woven coating	Fastening and combination options		CE 1 Items suitable for combination are not included in the scope of delivery. These are available to order separately from us.
CareFlow			Number of chambers		Reverse side	Twin holder	Elastic band/ Cuff	
UB 2000 US	2000	90/130		Tap		x		
UB 2000 USo/ SB 2000 USo	2000	90				x		
SB 1500 USo/ UB 1500 USo	1500	10/90				x		
UB 750 US	750	10/30/50	1	Tap	x	x	x	
UB 750 UF	750	50	1	Valve	x	x	x	
UB 500 US	500	10/30/50	3	Tap	x	x	x	

- Indication/Application Area:**
Intended for use on patients with a recumbent indwelling bladder catheter, per-cutaneous renal fistula or urinary sheath for collecting the drained urine by means of a urine drainage system.
- Draining the System:**
– If necessary, remove the drainage tap from the rear flap. Hold it over a suitable collecting container and open the bag by moving the non-return valve in the drainage tap. After the system has been completely emptied, push back the non-return valve, disinfect the drainage tap and insert it into the rear flap for safe and hygienic storage, if available.
- To empty the leg bag using the drainage valve, tilt the valve down 180°.
- Disposal:**
– The emptied bag system should be disposed of with the household waste in accordance with national regulations.
- Service Life/Duration of Use:**
– All urine drainage systems mentioned are intended for single use. If a system has been removed from its original packaging and connected to a catheter or urinary sheath, it must not be reconnected to the catheter (if the items have been separated). Closed urinary diversion systems should not be changed routinely at fixed intervals, but should be changed as necessary after being assessed by the physician or nursing staff, according to individual criteria.
- Packaging/Delivery:**
– Ensure that the packaging unit is unopened and undamaged and is not used for longer than indicated on the label.

- Warnings/Safety Precautions:**
Please note that urine drainage systems connected to an indwelling bladder catheter must not be pulled in any way in order to avoid disturbing the position of the catheter and impairing its drainage performance (RISK OF INJURY!)
- If the bag or medical device is defective and/or no urine enters the bag, discard it immediately and use a new one. If patients affected by conditions that require urinary diversion are fitted with our urine drainage systems, appropriate care/supervision must be ensured. The same applies to small children.
- Do not reuse, reprocess. Reuse, reprocessing may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in injury, illness or death of the patient. Reuse, reprocessing may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another.
- For device complaints, please contact your supplier. Please also report any serious incident that has occurred in relation to the device to the manufacturer and the competent authority at the user's place of residence.

Exclusion of Liability:
In general, warranty only applies if the connection of approved medical products from GHC to the product you have in hand is carried out correctly and in accordance with the instructions for use. If you prefer to connect these to products from other manufacturers, this is possible, but GHC assumes no liability and any warranty claims become void.

GHC assumes no liability whatsoever for injury or other damage resulting from use of this product other than in accordance with current GHC recommendations.

PLEASE NOTE; All information has been issued to the best of our knowledge at the time of printing. Subject to technical changes.

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