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GB CARE FLOW All Silicone Foley Catheter, Nelaton tip

DESCRIPTION

The product is made from silicone. It consists of shaft, drainage funnel, inflation funnel, straight tip (Nelaton tip) or Tiemann tip, irrigation funnel (if present), pre-loaded stylet (if present), balloon and valve. Spout and syringe are optional. The product can be used on perioperative patient, patient have urinary retention or bladder outlet obstruction, patient who need to measure the urine volume accurately.

FEATURES

- Made from 100% silicone
- Color coding for size identification
- With plastic valve
- DEHP-FREE

INTENDED USE / INDICATIONS

It is intended to be used by trained medical staff for drainage and/or irrigation of the patient bladder by inserting the catheter into the vesical cavity of the bladder through the urethra.

CONTRAINDICATIONS

- Acute urethritis
- Acute prostatitis
- Acute epididymitis

ADVERSE REACTIONS

- Reported adverse reactions associated with Foley Catheter are: Septicemia, Urethritis, Urinary tract infection and encrustations, bladder spasms, abdominal discomfort, etc.
- If difficulty is encountered aspirating the balloon with a syringe, a rare and infrequently reported event, and the leg of the catheter with the valve should be cut with a sharp scissors at the bifurcation or the balloon ruptured according to established procedures reported in medical literature. Should it be necessary to rupture the balloon, care must be taken to remove all fragments from the patient's bladder.
- Incorrectly positioned catheters can cause urethral damage if the balloon is inflated contrary to instructions within the urethra.
- Irritation of the urethral mucosa, blockage of the catheter due to encrustation and catheter induced infections are documented complications with some catheter materials and patients. The patient should be routinely monitored in accordance with accepted procedures and the catheter shall be removed after a suitable interval as determined by a trained medical staff or other suitably qualified personnel.
- Long term use of catheter has the risk of urinary tract infection, sphincter relaxation and urethral relaxation.

DESCRIPTION FOR USE

Preparation

- The following additional items are required for catheterization: Sterile field, Sterile gloves, Items required for cleaning the patient meatus (sterile, based on established techniques), Luer syringe with sterile water or sterile glycerin solution for balloon inflation (if not included), Empty syringe for balloon deflate(if not included), Sterile dressings
- ✓ Urine drainage device

- Before use, visually inspect the catheter to detect mechanical damages and check that no leakage occurs.
- Place male patients in supine position, female patients in lithotomy position.
- Wash and dry hands thoroughly.
- Using aseptic technique, remove the catheter from its pouch and place it on a sterile field.
- Put on sterile gloves and remove the catheter sleeve.
- Lubricate the tip and the shaft of the catheter with a water-based lubricant.
- Lean the opening of the urethra and the surrounding area, using established techniques.

Catheterization

- With the uncontaminated hand, pick up the catheter from the sterile field.
- Carefully insert catheter tip into the bladder (normally indicated by urine flow), and then a further 3 cm to ensure balloon is also inside the bladder.
- Inflate the catheter balloon with sterile water or sterile aqueous glycerin solution, at the discretion of the physician.
- Use the syringe or sterile glycerin solution to inflate the balloon with the required fluid volume (marked on the catheter funnel).
- Slowly retract the catheter until some resistance is felt to ensure that the balloon is correctly located within the bladder, before the bladder neck.
- Connect the catheter to a urine drainage device.
- Observe urine flow.

Catheter removal

- For removal of the catheter, deflate the balloon by inserting a Luer syringe into the valve. Release the syringe plunger and allow the balloon to deflate. Only use gentle aspiration, if necessary, to deflate the balloon. The balloon should spontaneously deflate.
- Discard the catheter according to hospital protocol.

Catheter Care

- Ensure local cleaning and hygiene protocols are followed to keep the catheter and meatus as clean as possible.
- Conduct hand hygiene immediately before and after any manipulation of the catheter and the urine drainage system. Wear disposable gloves when handling the system.
- Maintain unobstructed urine flow.
- Ensure that the catheter and collection tube do not kink.

LONG-TERM USING CARE

- The patient should wash his hands before and after any catheter operation.
- The area around the catheter is rinsed with a mild soap box at least twice a day.
- The patient taking 2 to 3 liters of water a day can help wash the catheter and keep the urine clear.
- The patient should regularly drain urine bags.
- The patient should always keep the drainage system closed.

WARNINGS / PRECAUTIONS

- Read all warnings and instructions before use. Improper use can result in serious or fatal illness or injury.
- Do not use if the sterile packaging is damaged or unintentionally opened before use.
- Do not use if it's damaged or irregularly shaped.
- Do not use after the expiry date.
- For single use only.
- Not compatible with MRI.
- Do not over inflate the balloon. Refer to outer unit pack or funnel of catheter for balloon capacity.
- Inflate catheter balloon only with sterile water or sterile glycerine solution.
- The pre-filled syringe with sterile water or sterile aqueous glycerine solution (if included) is for inflation of the balloon only. Not for injection.
- The empty syringe (if included) is for deflation of the balloon before the removal of the catheter.
- The product is intended for use by trained medical staff. Aseptic technique must be practiced.
- Do not use ointment and lubricants having petroleum base. Use water soluble lubricants.
- Do not clamp catheter shaft. It may damage catheter and prevent deflation.
- Do not use needle syringe to inflate balloon. Use a Luer syringe (Luer slip or Luer lock).
- Patients should be regularly monitored as determined by a physician.
- Ensure local cleaning and hygiene protocols are followed to keep the catheter and meatus as clean as possible.
- During Long-term use, the balloon inflation volume should be regularly monitored in case any clinical signs of deflation occur, such as bypassing of urine or urethral pain. If necessary, the balloon volume should be adjusted.
- The balloon must be completely deflated before catheter removal.
- The retention time of Silicone Foley Catheter should be decided by the medical care professionals according to the actual situation of using. The maximum indwelling time is 90 days. Suggest to use less than 30 days.
- Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in injury, illness or death of the patient. Reuse, reprocessing or resterilization may also create a risk

- of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.
- Precautions shall be taken when discarding device and disposal of the device shall be mad in accordance with applicable national regulations for biological hazardous waste.
- Dispose of product and packaging in accordance with hospital administrative and/or local government policy.
- Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the local competent authority of the user's place.

STORAGE CONDITIONS

- Store in a dark, cool and dry place.
- Protect product from moisture and excessive heat.
- Avoid prolonged exposure to ultraviolet, sunlight and fluorescent light.
- Store in manner preventing crushing.

CARE FLOW Transurethral SOFT Silicone Catheter, transparent

Article numbe	Type	Tip	Size	Color of Shaft	Balloon	Length
01/01/01/06/03/GHC	2-way	Nelaton	CH06	light red	03ml	31cm
01/01/01/08/05/GHC	2-way	Nelaton	CH08	black	05ml	31cm
01/01/01/10/05/GHC	2-way	Nelaton	CH10	grey	05ml	31cm
01/01/01/12/10/GHC	2-way	Nelaton	CH12	white	10ml	40cm
01/01/01/14/10/GHC	2-way	Nelaton	CH14	green	10ml	40cm
01/01/01/16/30/GHC	2-way	Nelaton	CH16	orange	30ml	40cm
01/01/01/18/30/GHC	2-way	Nelaton	CH18	red	30ml	40cm
01/01/01/20/30/GHC	2-way	Nelaton	CH20	yellow	30ml	40cm
01/01/01/22/30/GHC	2-way	Nelaton	CH22	violet	30ml	40cm
01/01/01/24/30/GHC	2-way	Nelaton	CH24	blue	30ml	40cm
01/01/01/26/30/GHC	2-way	Nelaton	CH26	pink	30ml	40cm

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