BRUGSINFORMATIONER DK

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
Rev. 003/08.2021

Hangzhou Primecare Medical Co., Ltd.
Room 408-409, Zancheng Center West
Shangcheng District
310008 Hangzhou
PEOPLE'S REPUBLIC OF CHINA

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INSTRUCTIONS FOR USE
GB

General Information:
 Dear user, please read all of the following information carefully before use. Please refer to the table for the specific product features of the individual systems.
Intended Purpose/Description:
 Urine drainage systems and bags for connection to sterile catheters. Non-sterile urine drainage systems and bags are only suitable for use with urinary sheaths and for connection to leg bags.

Device	Filling volume in ml	Tube length in cm	Bedside bag system	Leg bag system	Drainage tap/valve	Tube clamp	Pump ball	Non-woven coating	Twining holder	Cord	Twining holder	Elastic band	Support stocking/cuff	1 Items suitable for combination are not included in the scope of delivery. These are available to order separately from us.
CareFlow			closed	Number of chambers				Reverse side	Twin holder	Cord	Twin holder	Elastic band	Support stocking/cuff	
UB 4000 TUR	4000	120	x		Tap	x	x		x	x				
UB 4000	4000	120	x		Tap	x			x	x				
UB 2000	2000	120/90	x		Tap	x			x	x				
UB 2000A	2000	110			Tap	x			x					
UB 2000E	2000	90			Tap						x			
UB 1000 Day/Night	1000	30/60/90	x		Tap			x		x	x	x	x	
UB 750/750A	750	10/30/50		1	Tap			x		x	x	x	x	
UB 750F	750	50		1	Valve			x			x	x	x	
UB 500	500	10/30/50		3	Tap			x			x	x	x	

Indication/Scope of Application:
 Intended for use on patients with a recumbent indwelling bladder catheter, percutaneous renal fistula or urinary sheath for collecting the drained urine by means of a urine drainage system. The patient must be cared for in a way that is tailored to their individual needs. If the attending physician has given any previous instructions, then these must be observed.
Contraindication and Restriction of Use:
 We are not aware of any absolute contraindications. Restrictions of use may result from the clinical picture of the patient. Leg bag systems should not be used on patients with macrohematuria. Our urine drainage systems and their tubes must not come into contact with areas of the skin presenting pathological changes.
Application/Execution/Hygiene:
 Before manipulating the urine drainage systems in any way (connection or change of the system/urine drainage), please carry out hygienic hand disinfection (wash and disinfect hands). Please wear disposable gloves.
 – Take the bag out of the sterile packaging
 – Remove protective cap on the cone. If necessary, close non-return valve on bags at drainage tap or tilt valve (if open).
 – Separate the existing bag system from the catheter and disinfect the catheter funnel with suitable disinfectant. Observe the reaction time.
 – Connect the bag system to the catheter by inserting the cone into the catheter funnel.
 – Fasten the bedside bag system to the bed frame using twin holders and cord.
 – Leg bag: For optimal drainage, ensure that the bag is placed below the bladder.
 – Please take care that the tubes of our urine drainage systems, especially the leg bags, are not kinked to avoid accidentally disrupting the drainage.

Disposal:
 – The emptied bag system should be disposed of with the household waste in accordance with national regulations.
Service Life/Duration of Use:
 – All urine drainage systems and bags mentioned are intended for single use. If a system has been removed from its original packaging and connected to a catheter or urinary sheath, it must not be reconnected to the catheter if the items have been separated from one another. Closed urinary diversion systems should not be changed routinely at fixed intervals, but should be changed as necessary after being assessed by the physician or nursing staff, according to individual criteria.
Packaging/Delivery:
 – Each sterile item is delivered in individual packaging. Sterility is only guaranteed as long as the packaging unit is unopened and undamaged, but not longer than is indicated on the label.

Warnings/Precautions:
 Please note that urine drainage systems and bags connected to an indwelling bladder catheter must not be pulled in any way in order to avoid disturbing the position of the catheter and impairing its drainage performance (RISK OF INJURY)! If the bag or medical device is defective and/or no urine enters the bag, discard it immediately and use a new one. If patients affected by conditions that require urinary diversion are fitted with our urine drainage systems or bags, appropriate care/supervision must be ensured. The same applies to small children.
Exclusion of Liability:
 In general, warranty only applies if the connection of approved medical products from GHC to the device you have in hand is carried out correctly and in accordance with the instructions for use. If you prefer to connect these to devices from other manufacturers, this is possible, but GHC assumes no liability and any warranty claims become void.
 GHC assumes no liability whatsoever for injury or other damage resulting from use of this device other than that which is in accordance with current GHC recommendations.
 PLEASE NOTE: All information has been issued to the best of our knowledge at the time of printing. Subject to technical changes.

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Manufacturer	Importer	Catalogue number	Do not resterilize	Keep away from sunlight	Do not re-use
Authorized representative in the European community	Use-by date	Medical device	Single sterile barrier system	Keep dry	Consult instructions for use
Date of manufacture	Batch code	Sterilized using ethylene oxide	Do not use if package is damaged	Lower limit of temperature	Contains or presence of phthalates

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